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Manual instruction Blood Pressure Monitor

MM – 210 ESATTO

1.	Before Using This Device	25
2.	Important Safety Information and Precautions	25
3.	Specifications.	26
4.	About Blood Pressure.....	26
5.	About blood pressure monitor	27
6.	Buttons.....	27
7.	Display.....	27
8.	Battery installation	28
9.	Power supply from USB port	28
10.	Putting on the cuff	28
11.	Body positioning during measurement.....	29
12.	Measurements notes	29
13.	Setting of user.....	29
14.	Turn on the device	29
15.	Setting of Year, Date, Time, 24-hour system/ 12-hour system, pressure unit	29
16.	Measurement	30
17.	Function for the detection of cardiac disorders - arrhythmias.....	31
18.	Using Memory-Recall and Memory-deleting	31
19.	Special function - determining the average of the last 3 measurements	31
20.	Memory view	31
21.	Remove measurement results from the device's memory	32
22.	Troubleshooting.....	32
	Terms of warranty	32
	Dziennik pomiarów / Measurement log	34

Thank you for purchasing our product. Please read manual instruction before use.

1. Before Using This Device

INDICATIONS FOR USE: This blood pressure monitor is intended to be used to measure blood pressure and pulse rate from upper arm by using the oscillometric method at home. This device is intended for use on adults only and is not for use on groups such as neonatal babies, children, pregnant women, or pre-eclamptic patients. It can't be used while the arm has bleeding or wound to avoid the blood flowing from the wound in pressurizing.

2. Important Safety Information and Precautions

MEDICAL DISCLAIMER: This device and manual are not meant to be a substitute for advice provided by doctors or other medical professionals. Contact your physician for interpretation of measurements, or if you have or suspect you have a medical issue.


- Contact your physician if you suspect you have a medical issue.
- Not for use by children, pregnant women, or pre-eclamptic patients.
- If you suffer from arrhythmia, diabetes, blood circulation problems, or apoplexy, only use this device under a physician's care.
- Prolonged over-inflation of the monitor will result in harmful injury to the patient.
- Do not use this device over a wound, as this can cause further injury.
- Too frequent measurements can cause injury due to blood flow interference.
- This device possible cause improper result due to external interference, such as acceleration during transport or transport in general.
- Do not attempt to modify this device in any way.
- Store this device in a cool, dry place. Do not subject this device to extreme temperatures, humidity, or sunlight. This device might not meet performance specifications if stored or used outside the ranges specified in the "Specifications" section. Keep away from rain.
- Do not immerse this device in water or clean with cleaning products, alcohol, or solvents. Carefully follow cleaning instructions provided.
- Remove batteries if this device will not be used for three or more months.
- Ensure that the device is used in the environment specified in the declaration of conformity in this manual, otherwise, may result in improper operation.
- Use of this device adjacent to or stacked with other device should be avoided because it could result in improper operation.
- This electrical medical equipment requires specific precautions regarding electromagnetic compatibility. It must be installed and used according to the electromagnetic information.
- Do not use this device where flammable gases or liquids are present.
- Never drop the device

Technical description is contained in this instruction manual.

3. Specifications.

Model	MM 210 Esatto/BP101k
Measurement Method	Oscillometric
Display	LCD
Measurement Range	Pressure: 0 to 299 mmHg Pulse: 40 to 180 beats/minute
Accuracy	Pressure: Within 3 mmHg Pulse Rate: Within 5% of reading
Inflation	Automatic inflation by pump
Deflation	Automatic rapid deflation
Pressure Detection	Semiconductor pressure sensor
Memory	persons X 90 measurements each
Power supply	3 AA alkaline batteries (not included)
Battery	About 6 months when used once a day
Operating Environment	Temperature 50°F (10°C) to 104°F (40°C) Humidity 15% to 90% (noncondensing) Barometric ure 105 kPa to 80 kPa (790 mmHg to 600)
Transport Storage Environment	Temperature -4 °F (-20°C) to 131°F (55°C) Humidity 0% to 95% (noncondensing)
Weight of Main Unit	115 g (without batteries)
External Dimensions	103 mm (H) X 77,5mm (W)X 40,8 mm
Measurable Circumference of arm	22~36cm
Electric Shock Protection	Internal power supply appliance type BF
Contents	Main unit, cuff, manual

4. About Blood Pressure

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/ stethoscope auscultatory method, within the limits prescribed by European Society of Hypertension International Protocol revision 2010 for the validation of blood pressure measuring devices in adults. An irregular heartbeat (IHB) rhythm is defined as a rhythm that is more than 25% slower or 25% faster from the average rhythm detected while the device is measuring the systolic blood pressure and the diastolic blood pressure. When the device detects an irregular rhythm two or more times during the measurement, the Irregular Heartbeat Symbol will be displayed with the measurement values.

5. About blood pressure monitor


6. Buttons

„Start/Stop” button - Universal button, start of measurement.

„Memory button” - it is used to view memory, change the settings parameters.

„Set button” - Switch between users, moving between settings.

7. Display

8. Battery installation

- a. Take off the battery cover on the back of this device.(Figure 1)

Figure 1

- b. Insert three AA alkaline batteries according to the + and – markings. (Figure 2)

Figure 2

- c. Snap the battery cover on, ensure it clicks into place.

9. Power supply from USB port

Plug the micro USB cable into the micro USB port on the device as shown (Figure 3).

Figure 3

10. Putting on the cuff

The cuff should be put on the left hand. Distance from the cuff to elbow bend is 1-2 cm.

11.Body positioning during measurement

Before measuring blood pressure, you should sit comfortably straighten and relax the body.

It is necessary to turn the hand with the inner side up and place the hand at the same height as the heart.

12.Measurements notes

- After measurements, the blood pressure monitor turns off after 2 minutes after pressing the last button.
- To stop the measurement, press the START/STOP button.
- Between measurements, keep a minimum break of 5 minutes and try to maintain one body position and breathe evenly and calmly.
- Carefully attach the cuff.
- Do not tilt or tighten the muscles of the hand and hand during measurements.
- Try to maintain one body position.
- Do not talk or move during measurements.
- The result will be automatically stored in the previously selected memory when the measurement is complete.

13.Setting of user

With device turned off, press and release “SET” button to show the current user.Press “SET” button to toggle between User 1 and User 2. To select user, press “START/STOP” to confirm and then the device turns off.

14.Turn on the device

1. Press Start/Stop.
2. When you press the START/STOP button, the device is ready.
3. The blood pressure monitor has an automatic shut-off function and will turn off 2 minutes after the last press of the button.

15.Setting of Year, Date, Time, 24-hour system/ 12-hour system, pressure unit

- a. With device turned off, press and hold “SET” button until flashing year digits occur.
- b. Press “MEM” to increase one year, or press and hold “MEM” to increase years. Press “SET” to confirm and then flashing month digits occur.
- c. Press “MEM” to increase one month, or press and hold “MEM” to increase months. Press “SET” to confirm and then flashing day digits occur.

- d. Press “MEM” to increase one day or press and hold “MEM” to increase days.
Press “SET” to confirm and then flashing hour digits occur.
- e. Press “MEM” to increase one hour or press and hold “MEM” to increase hours.
Press “SET” to confirm and then flashing minute digits occur.
- f. Press “MEM” to increase one minute or press and hold “MEM” to increase
minutes.
Press “SET” to confirm and then either “12” digits occur for 12-hour system or
“24” digits occur for 24-hour system.
- g. Press “MEM” to toggle between “12” for 12-hour system and “24” for 24-hour
system. Press “SET” to confirm and then pressure unit digits occur.
- h. Press “MEM” to toggle between “mmHg” and “kPa”. Press “SET” to confirm
and flashing year occur again, or press “START/STOP” to confirm and then the
device turns off.
- i. Setting is completed.

Note: Pressing “START/STOP” will confirm the setting and then the device turns off at any time.

16. Measurement

- a. This device should be used on bare skin. If necessary, loosely roll up your sleeve
to expose arm. Note that your arm cannot be tightened by your sleeve.
- b. With the hose down the inside of your arm and the bottom of the cuff approx.
1~2cm above your elbow, wrap the cuff around your arm. Secure the cuff with a
finger space between the cuff and your arm by attaching hook-and-loop fastener.
- c. Comfortably seated
 - Legs uncrossed
 - Feet flat on the floor
 - Back and arm supported
 - Relax as much as possible and do not talk during the measurement.
- d. While sitting upright in a chair with your feet flat on the floor, press
“START/STOP”.

- Symbols on screen will flash and cuff will begin to inflate. Do not move, activate arm muscles, or talk during measurement.
- e. Measurement can be stopped at anytime by pressing “START/STOP”.
 - f. When a pulse is detected, a heart symbol will flash on screen.
 - g. When measurement is complete, the cuff deflates automatically and your systolic pressure, diastolic pressure, and pulse rate will display on screen.
 - h. Turn off this device by pressing “START/STOP”, remove the cuff from your arm by detaching hook-and-loop fasteners. If this device is accidentally left on, it will turn itself off after two minutes.

NOTE: This device will remain set on the same user after you turn it off. To change to the other user, turn this device off, then press and release “SET” to access user setting. Either the User 1 or User 2 icon will appear. Press “SET” to switch to the other user icon and press “START/STOP” to confirm the user and then the device turns off.

17. Function for the detection of cardiac disorders - arrhythmias

The device has the function of detecting cardiac arrhythmias. After performing the measurement, an icon will appear on the display. Speaking, moving affects the measurement result which is why the device can read it as not a regular heart rhythm. If this icon appears, we recommend repeating the test. However, if this icon appears frequently, you should consult a doctor.

18. Using Memory-Recall and Memory-deleting

This device will store up to 90 measurements each for two users. When more than 90 measurements are taken, the oldest measurement will be replaced.

- To access Memory-Recall, press and release “MEM” while this device is off.
- First, it will display the average of your three most recent measurements, indicated on screen as M AVG (Measurement Average).
- Press “MEM” again to see the most recent measurement. Press “MEM” repeatedly to scroll towards the oldest measurement. On screen, M01 indicates the most recent measurement, with higher M numbers indicating older measurements.
- If you desire to delete all stored measurements, press and hold “SET” under the Memory-recall status until “dEL no” is displayed with “no” flashing. Press and release “SET” to switch to “dEL YES” with “YES” flashing. Press “START/STOP” with “dEL YES” to confirm clearing the memory and then blank entry appears on screen, or press “START/STOP” with “dEL no” to go back to the Memory-recall status without clearing the memory.
- Press “START/STOP” to exit.

19. Special function - determining the average of the last 3 measurements

After each measurement, the result is stored in the device's memory. When you press the MEMORY button once, the average of the last 3 measurements is shown on the display. This feature is very useful for regular measurements at the same time of day, allowing you to even more accurately illustrate your health.

20. Memory view

To view measurements stored in memory:

1. Press the MEMORY button, the display will show the average of the last 3 measurements.
2. By pressing the MEMORY button, you can view subsequent measurements.

21. Remove measurement results from the device's memory

Press “START/STOP” with “dEL YES” to confirm clearing the memory and then blank entry appears on screen, or press “START/STOP” with “dEL no” to go back to the Memory-recall status without clearing the memory.
Press “START/STOP” to exit.

22. Troubleshooting

If you see any of these symbols displayed on screen, perform the corresponding correction:

Err Error during inflation

Correction: Adjust cuff, ensuring it is wrapped properly, and repeat measurement.
Do not move or talk during measurement. Ensure cuff is level with heart.

Lo Low battery

Correction: Replace batteries following instructions in the Battery installation section

If you have the following issues without an error code, perform the corresponding correction:

No Power Low batteries or incorrect battery placement *Correction:* Ensure batteries were installed according to the + and - shown on this device. If they were, batteries need replacement.

Measurement Too High / Too Low Possible incorrect use of this device

Correction: Adjust cuff, ensuring it is wrapped properly, and repeat measurement.