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czynne w godz. 7.30 -15.30

1.Feature

- Adopts reliable sensor made in German.
- All-new & patented probe design ensures high accuracy.
- Excellent adaption to ambient temperature. Accurate and reliable even under complicated surrounding.
- Body mode and object mode available.One switch to change mode.
- Beeper for fever or high temperature (alarm value can be defined by user).
- Automatic recall of last reading.
- Large LCD with back-lit.
- C/F reading available.
- Automatic shut off.

2.Safety precautions

- Using this thermometer is not intended as a substitute for consultation with your physician. It is dangerous for users to perform a self-evaluation and self-treatment based on the measuring result. please follow doctor' s instruction.
- Keep the thermometer out of reach of children ,please consult the doctor at once if the children accidental swallow of battery or other component .
- Don't throw the battery into fire.
- The device is precision instrument, don' t drop, tramp or impose any vibration or impact on the thermometer.
- Do not touch the lens of the probe with your fingers and disassemble the device by yourself.
- Please make sure your forehead is clean before measuring forehead temperature.
- Please stay still indoors about 30 minutes after you do some exercise, like eating and bathing before measuring.
- Please place the thermometer indoors for about 30 minutes if ambient temperature varies a lot before using.
- When measure somebody continuously, the temperature should be measured every minute, if you need to measure yourself continuously for a short time, It is a normal phenomenon that there are some slight errors when you read the temperature. We recommend that you measure yourself continuously maximum of three in a unit of time then figure out the average and choose it because the temperature of the human will conduct to the thermometer, it may affect the accuracy of measurement.
- please collect the record of Individual temperature under the good condition of body in usual days as a reference for checking fever or not.
- Do not measure the sites of scarred tissue or tissue compromised by skin disorders because they will affect the accuracy of measurement.
- Do not measure the site of forehead temperature if patients has trauma on forehead.

- Do not measure if patient is treated with certain drug therapies because body temperature may increase in the drug within the effort time limit.
- Do not immerse the device into water or any other liquid, and not directly sunlight exposure.
- Do not use a mobile or cordless phone near the thermometer when measuring.
- Please don't measure body temperature in strong electromagnetic interference environment (such as microwave, high frequency equipment operation environment) to ensure the accuracy of measurement data.
- Do not disassemble, repair, or modify the unit.
- This thermometer only a personal device, please do not share with others.
- Do not touch the battery output when measuring.
- Please store the thermometer according to the technical specifications.
- The materials (ABS) of external contact with patient has passed the ISO 10993-5 and ISO 10993-10 standard test, no toxicity, allergy and irritation reaction. They are compliant with the MDD requirements based on the current science and technology, and other potential allergic reactions are unknown.
- The patient is an intended operator. The patient can measure, read data and replace battery under normal circumstances and maintain the device and its accessories according to the user manual
- Don't use this thermometer for other purposes.
- It is forbidden to leave the product exposed to any chemical solvent, direct sunshine or high temperature in case of damaging the product or the battery.
- Do not measure while talking on the phone.
- Please report to MANUFACTURER if any unexpected operation or events occurs.

3. Battery replacement


1. Open the battery lid and take out exhausted battery.
2. Put into 2 AA alkaline batteries and close up battery lid. After new battery is installed, thermometer will give "Beep.Beep". If there's no beeps, check if the positive and negative pole is correct (see fig.10.1).

Notice

- 1). Take out battery in case the thermometer is not used for long period.
- 2). Dispose battery according to local regulations. Don't put the battery to fire

NOTICE:


- Please observe the related national laws of disposing the abandoned battery and don't litter to the garbage can.
- Please take out the battery if the device is not used for long periods of time.
- Please don't put the battery into the fire.
- To protect the environment, dispose of empty battery at your retail store or at appropriate collection sites according to national or local regulations.


4.Specification

Model:	HW-3/MM007	
Measurement method	Non-contact	
Effective distance	5CM ~ 8CM	
Range	Human body	35.5°C ~ 42.9°C (95.9°F ~ 109.2°F)
	Object	0°C ~ 118°C (32°F ~ 244°F)
Accuracy	Human body	±0.2°C/0.4°F
	Object	±1.0°C/1.8°F
Resolution	0.1°C/°F	
Working condition	10°C ~ 40°C (50.0°F ~ 104°F) RH ≤ 95% Non-condensing.	
Storage condition	-25°C ~ 55°C (-13°F ~ 131°F) RH ≤ 95% Non-condensing.	
Power supply	d.c. 3V 2 AA alkaline Battery	
Power consumption	When off ≤ 10uW	
	When measurement ≤ 30mW	
Power level indicator	Indication for low power level.	
Memory	Automatic recall of last reading.	
Display	LCD back-lit display	
Reading scale	Celsius or Fahrenheit	
Automatic shut off	In 30 seconds	
Dimensions	93mm × 153mm × 41mm	
Weight	125g	

5.Illustration

6.Display and icons

- 1.Celsius reading.
- 2.Fahrenheit reading.
- 3.Memory - Temperature value of last measurement.
- 4.Object mode.
- 5.Body mode.
- 6.Result
- 7.Beeper
- 8.Battery level

7.Maintenance & tips

•Make sure the sensor and probe cavity is clean otherwise it will affect accuracy.
Cleaning method for probe:

- 1.Use the cotton stick or soft cloth with water or alcohol to wipe the casing.
 - 2.Use the cotton stick or soft cloth with alcohol to wipe the sensor surface or probe cavity gently. Don't use thermometer before alcohol is vaporized.
- Read this manual book thoroughly before use. Make sure battery is well installed.
 - It is not allowed to put the thermometer in any liquid or expose to strong sunlight or extremely low temperature..
 - Strong crash or hit to the product will cause its damage. Do not dismantle this thermometer by yourself.
 - Keep the thermometer from children' s reach.
 - Do not use the thermometer under circumstance of strong electromagnetic interfere.
 - The measurement results are probably fluctuating due to improper measurement ways. Please practice adequate measurements in order to improve your skill.
 - The measurement results can not supersede a doctor' s diagnosis.
 - Special maintenance is unnecessary for this thermometer. Please contact distributor or manufacturer in case of malfunction.

Care and cleaning

- The probe tip and lens are the most delicate part of the thermometer. It has to be clean and intact to ensure accurate readings.
- Wipe the surface gently with a cotton swab or soft cloth moistened with alcohol until the alcohol completely dried out.
- If the lens is damaged, contact the distributor.

Clean the unit body:

- Use a soft, dry cloth to clean the thermometer display and unit body.
- If very dirty, use a soft cloth with alcohol to cleaning.

NOTES:

- Don not use abrasive cleaners.
- Don' t use other non-recommended methods to perform disinfection.
- Non-waterproof, don' t use the abrasive cleaner to clean the product or drop the thermometer into water or other liquid.

Maintenance

- 1) We do not authorize any institution or individual to maintain and repair the product. If you suspect that the products have any questions, please contact the manufacturer or distributor to handle the case.
- 2) The user must not attempt any repairs to the device or any of its accessories. Please contact the retailer for repair.
- 3) Dismantling equipment by unauthorized agencies is not allowed and will terminate any claim to warranty.


Calibration

The thermometer has been initially calibrated at the manufacturing process. If this thermometer is used according to the use instruction, periodic re-adjustment is not needed. If you have question about the accuracy of measurement at any time, please contact distributor or manufacturer, the contact information see last page.

Storage

- 1) Don't put the thermometer under the sunshine, high temperature and moist environment or someplace which maybe get in touch with fire or is vulnerable to vibration.
- 2) Take out the battery if don't need to use the device in a long time.

8.Operation Instruction

- 1.Body/Object - To switch measurement mode between human body and object
- 2.Menu- Activates menu and save setting
- 3.,+ " - Increases parameter or select relevant parameter when setting parameter
- 3.,- " - Reduces parameter or select relevant parameter when setting parameter

9.Preparation

Check battery

- 1.Replace the batteries to ensure power supply if the thermometer's LCD display low voltage icon.

2.Check sensor

If the sensor have pollution and spray, please clean it. (The cleaning method see the chapter 8 Care and Cleaning for detailed.)

If the sensor lens is damaged, please stop using.

- 3.Please put the thermometer into the measurement environment(16°C~40°C) for 30 minutes to gain more accurate measuring result.

- 4.Accuracy of unexpected fluctuations in ambient temperature may decrease the

measurement results. When the thermometer at the same measurement position display different ambient temperature, or test temperature in front of the air conditioner, it will not be able to obtain accurate results.

5.If you want to measure forehead temperature, please make sure the forehead is naked and clean to ensure the accuracy of measurement.

10.Settings

User can change reading scale between Celsius or Fahrenheit, switch on or switch off beeper; set alarm value, set deviation value and change measurement mode between human body or object. Please refer to below chart.

Measurement mode setting:

When thermometer is on, it displays current measurement mode (fig 7.1).

Switch the "Body/Object" button to change measurement mode (Fig 7.2).

Fig 7.1


Fig 7.2

Reading scale setting: F-1

Activates setting menu by pressing "MENU" button. "F-1" will be displayed. Press "+" to select Fahrenheit "°F" as temp unit (icon "°F" is twinkling, Fig 7.3); or press "-" to select centigrade "°C" as temp unit (icon "°C" is twinkling, fig 7.4).

Fig 7.3

Fig 7.4

Beeper setting: F-2

Press "MENU" again to save the setup and "F-2" will be displayed (Fig 7.5). Beeper is defaulted as "OPEN". Press "+" to switch on beeper and "OPEN" will be displayed (Fig 7.6). Press "-" to switch off beeper and "CLOS" will be displayed (Fig 7.7).

Fig 7.5

Fig 7.6

Fig 7.7

Alarm value setting: F-3

Press MENU again to save the setup and "F-3" will be displayed (Fig. 7.8). Defaulted value "38.0°C" will be displayed (Fig. 7.9). Press "+" or "-" to increase or decrease alarm value (Fig. 7.10). Value will be increased or decreased by 0.1°C step by step. Fig. 7.8 Fig. 7.9 Fig. 7.10

Fig 7.8

Fig 7.9

Fig 7.10

Deviation value setting: F-4

Press MENU again to save the setup and "F-4" will be displayed (Fig. 7.11). Defaulted deviation value is 0.8°C (Fig 7.12). When defaulted deviation value is twinkling, press "+" or "-" to increase or decrease deviation value (Fig. 7.13). Deviation value will be increased or decreased by 0.1°C step by step. Fig. 7.11 Fig. 7.12 Fig. 7.13

Fig 7.11

Fig 7.12

Fig 7.13

Save and exit (automatic shut off):

Press MENU to save the setup and "SAVE" will be displayed (Fig. 7.14). All setup will be saved and thermometer will automatically shut off (Fig. 7.15).

Fig 7.14

Fig 7.15

After activating setup menu, users may select specific functions (F1/F2/F3/F4) directly to do specific setting.

New setup will not be saved if the thermometer is shut off unexpectedly during the procedure of setting.

After activating setup menu, the thermometer is unable to measure temperature. There's no response when pressing measurement button.

Notice:

1. Temperature under human body mode is obtained from dynamic compensation of environmental temp and forehead surface temp.

2.Object temperature mode is to test surface temperature of an object. The temperature get from forehead under this mode id merely temperature of forehead surface but not body temperature.

3.Deviation setting is able to adjust measurement value from 0°C to 3.0°C according to testing distance, surrounding temperature, skin difference, etc. Defaulted value is 0.8°C.

For example: if the tested body temperature is 36.2 while the temperature is supposed to be 37.0°C, "F-4" to set up the value and adjust to 1.6°C, then the measurement value will be 37.0°C.

11.Measurment

Body temperature Measurment

- Press measurement button to turn on thermometer and it displays boot screen (Fig.7.16). After POST and two beeps, it will display value of last reading and be ready for measurement (Fig. 7.17).
- Make sure the thermometer is under body mode.
- Keep distance at 5cm to 8cm from upper eyebrows to the probe (Fig, 7.18). Press measurement button and when it gives a "beep" measurement is finished and value will be displayed (Fig. 7.19). If measurement value is exceeding alarm value (Defaulted value is 38°C), it gives "beep. beep. Beep" as an indication.
- After measurement, if the thermometer is idle in 30 seconds, it will display "OFF" (Fig. 7.20) and gives a "beep" and shut off automatically


Fig 7.16


Fig 7.17


Fig 7.19

Fig 7.20


Fig 7.18

Notice:

- 1.Keep the sensor and probe cavity clean before and after use.

- 2.To ensure the accuracy of measurement, it is recommended to start measurement after ten minutes when carrying the thermometer to a new environment.
- 3.Wait for 10 minutes to measurement body temperature after measuring extremely high temperature or extremely low temperature objects.
- 4.Wait for 5 minutes to start a measurement when a measuring target (object or human) is from an environment with enormous difference in temperature.
- 5.Breeze, water, sweating, cosmetic on forehead may affect measurement. Do not measure body temperature in 30 minutes after exercise, bath or meals.

Object temperature

- Press measurement button to turn on thermometer (Fig 7.21).
- Make sure the thermometer is under object mode.
- Keep vertical distance at 5cm to 8cm from object to measurement probe. Press measurement button and when it gives a “beep” measurement is finished and value will be displayed (Fig 7.22).
- After measurement, if the thermometer is idle in 30 seconds, it will display “OFF” (Fig 7.23) and gives a “Beep” and shut off automatically.


Fig 7.21


Fig 7.22

Fig 7.23

No response when pressing measurement button

1. Take out and reassemble battery.
2. Check if the thermometer is under menu setup. In procedure of menu setting, thermometer is unable to measure and therefore no response.

No display or improper display

1. Take out battery and install battery again.

No beeper

1. Check if the beeper is switched off.

Shut off right after switching on

1. Check battery level or take out and install the battery again.

12.Normalized symbols

	Follow operating instructions
	BF type applied part
	Caution

	Disposal in accordance with Directive 2012/19/EU (WEEE).
 1639	Complies with the European Medical Device Directive (93/42/EEC and amended Directive 2007/47/EC. Notified Body is SGS.
	Manufacturer information
	Authorized representative in the European Community
IP22	IP code of the device: this device' s grade of against ingress of solid foreign objects -- $\geq 12.5\text{mm}$ diameter (and the against access to hazardous parts with finger); the grade of waterproof is dripping (15° tilted).
	Batch code
	Date of manufacture
	Medical device

1) The digital thermometer needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.

2) Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can effect this equipment and should be kept at least a distance $d=3.3\text{m}$ away from the equipment. (Note: As indicated in Table 6 of IEC 60601-1-2:2007 of ME EQUIPMENT, a typical cell phone with a maximum output power of 2W yields $d=3.3\text{m}$ at an IMMUNITY LEVEL of 3V/m).

Mescomp Technologies S.A

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03-948 Warszawa, mescomp.pl

Tel. +48 (22) 636 20 00 mesmed.pl

Service Office

tel. +48 (22)838 63 38

Working hours 7.30-15.30